



Phoenix Health s.r.o

THE PROCESS OF PRODUCING MIDSTREAM PREGNACY RAPID TEST BY PHOENIX HEALTH S.R.O

Uncut sheet is ready to be processed.

1 - Cutting

Laminated plastic cards are cut into the size of 4.8mm x 95mm with a cutter machine.



2- Pocketing

2-1 Cut test strips are then assembled into the complete plastic Cassette by workers.



📍 Ing. Štefana Višňovského 211/6
Banská Štiavnica 969 01
Slovakia

☎ +421 (907) 773222



Phoenix Health s.r.o

2-2 Assembled cassettes will be put in a aluminum pouches with a silica gel pocket.



2-3 Quality control is implemented by direct supervision of operators.

3- Sealing and Printing

Individual Cassettes are sealed together with silica gel using aluminum foil pouch. The pouch is clearly labeled with lot number and an expiry date.





Phoenix Health s.r.o

4. Packaging Process

4- 1 Sealed and printed aluminum pouches containing the cassette kit are then packed into boxes.



4- 2 Final kit inspection of the outgoing kits is carried out by QC technicians prior to shipment.

5- In-process Control of Production Process

The in-process control steps are introduced at strategic steps to monitor the quality of products. The in-process controls performed within the various levels of production, are in accordance with in-house standard operating procedure.

THE SPACES FOR THIS PRODUCTION LINE DEPENDS ON THE QUANTITY THAT IS PLANNED TO PRODUCE. FOR STARTING AND SMALL SCALE 80 SQUARE METTER IS NECESSARY.

📍 Ing. Štefana Višňovského 211/6
Banská Štiavnica 969 01
Slovakia

☎️ +421 (907) 773222